

Installation


WITH NO PHONE

If you should experience any problems when first setting up your unit or if you have any questions, please call us at: **1-800-906-0872**

Telephone Jack



Emergency Unit



Step 1 Remove unit from box.

2 Place unit near an existing wall jack.

3 Using the provided phone cord, plug the first end into the existing wall jack. Plug the second end of the cord into the unit's adaptor, in the port labeled **"wall jack (A)"**.



4 Plug the transformer into an open outlet not controlled by a wall switch.



5 Turn unit on it's side. Locate the **"on/off"** switch, which can be found in the hole under the phone cords. The **"on/off"** lever has a white dot of paint on the end. When you first receive the unit the switch is set at the **"off"** position. Slide it to the **"on"** position.



6 Call **1-800-325-1787** to put your account on **TEST MODE**. After you hang up with the operator, press your emergency transmitter button. The unit will say **"EMERGENCY"** approximately 5-6 times and then say **"EMERGENCY REPORTED"**. The operator will come on after approximately 15-30 seconds and let you know that the signal was received. Talk back without picking up the phone. The operator will verify your address and answer your immediate questions.

7 In the event of an actual medical emergency, simply push the medical button. The operator will come over the speaker just as they did when you tested the system.

8 Congratulations! Your system is now on-line and ready to use.

