


# Traditional Installation

## USING BASIC LANDLINE

If you should experience any problems when first setting up your unit or if you have any questions, please call us at: **1-800-906-0872**


Telephone Jack



Emergency Unit



Telephone



**Step 1** Remove unit from box.

**2** Place unit near an existing wall jack and phone.

**3** Locate the cord that stretches from the wall jack to the phone. Unplug the end that goes into the phone and plug it into the unit's adaptor, in the port labeled "**wall jack (A)**".



**4** Using the provided phone cord, connect your existing phone into the unit's adaptor, in the port labeled "**phone (B)**".

**5** Plug the transformer into an open outlet not controlled by a wall switch.



**6** Turn unit on its side. Locate the "**on/off**" switch, which can be found in the hole under the phone cords. The "**on/off**" lever has a white dot of paint on the end. When you first receive the unit the switch is set at the "**off**" position. Slide it to the "**on**" position.



**7** Call **1-800-325-1787** to put your account on **TEST MODE**. After you hang up with the operator, press your emergency transmitter button. The unit will say "**EMERGENCY**" approximately 5-6 times and then say "**EMERGENCY REPORTED**". The operator will come on after approximately 15-30 seconds and let you know that the signal was received. Talk back without picking up the phone. The operator will verify your address and answer your immediate questions.

**8** In the event of an actual medical emergency, simply push the medical button. The operator will come over the speaker just as they did when you tested the system.

**9** Congratulations! Your system is now on-line and ready to use.

## CORDLESS PHONE

Cordless phone users must plug into main base that's plugged into unit, not a satellite base.

